



Sustainable Solutions for health, productivity and the environment.

I T E M 1 4 7 3 2 8 2

# ENVIRO CARE® NEUTRAL DISINFECTANT EFFICACY DATA

## EPA REG #47371 - 131 - 527

Enviro Care Neutral Disinfectant is bactericidal against several gram-negative pathogens according to the AOAC Use Dilution Test (Current Edition), modified in the presence of 5% organic serum and 400 ppm water hardness (calculated as CaCO3). Dilution 1:64 (2 oz. per gallon of water).

- |                         |                                              |
|-------------------------|----------------------------------------------|
| Pseudomonas aeruginosa  | Klebsiella pneumoniae                        |
| Proteus vulgaris        | Fusobacterium necrophorum                    |
| Salmonella choleraesuis | Enterococcus faecalis (vancomycin resistant) |
| Proteus mirabilis       | Acinetobacter calcoaceticus                  |
| Enterobacter aerogenes  | Bordetella bronchiseptica                    |
| Serratia marcescens     | Chlamydia psittaci                           |
| Enterobacter cloacae    | Legionella pneumophila                       |
| Shigella flexneri       | Listeria monocytogenes                       |
| Salmonella typhimurium  | Pasteurella multocida                        |
| Shigella sonnei         | Salmonella enteritidis                       |
| Salmonella typhi        |                                              |
| Escherichia coli        |                                              |

Enviro Care Neutral Disinfectant is bactericidal against several gram-positive pathogens according to the AOAC Use Dilution Test (Current Edition), modified in the presence of 5% organic serum and the 400 ppm water hardness (calculated as CaCO3). Dilution 1:64 (2 oz. per gallon of water).

- Community Associated Methicillin-Resistant Staphylococcus aureus - (CA-MRSA) (NRS384) (USA300)
- Community Associated Methicillin-Resistant Staphylococcus aureus - (CA-MRSA) (NRS123) (USA400)
- Staphylococcus aureus-Methicillin resistant (MRSA)
- Staphylococcus epidermidis
- Streptococcus pyogenes
- Staphylococcus aureus-Vancomycin Intermediate Resistant (VISA)
- Streptococcus faecalis

Enviro Care Neutral Disinfectant is virucidal against HIV-1 (Aids) as tested according to the official HIV-1 Virucidal Qualification Test, modified in the presence of 5% organic serum and 400 ppm water hardness (calculated as CaCO3). Dilution 1:64 (2 oz. per gallon of water).

Enviro Care Neutral Disinfectant is virucidal against the following viruses as tested according to the virucidal Qualification Test, modified in the presence of 5% organic serum and 400 ppm water hardness (calculated as CaCO3). Dilution 1:64 (2 oz. per gallon of water).

- |                                         |                                                                          |
|-----------------------------------------|--------------------------------------------------------------------------|
| Avian polyomavirus                      | Aedonovirus, Type 4                                                      |
| Avian Influenza (H5N1)                  | Respiratory Syncytial Virus (RSV)                                        |
| Hepatitis-B Virus                       | Human coronavirus                                                        |
| Feline Leukemia                         | Newcastle Disease                                                        |
| Hepatitis-C Virus                       | Rubella (German Measles)                                                 |
| Canine Distemper                        | Rotavirus [WA]                                                           |
| Herpes Simlex Type II                   | Influenza Virus (in association with the organism Influenza A/Hong Kong) |
| Rabies                                  | SARS Associate                                                           |
| Feline Picornavirus                     | Coronavirus (SARS) (cause of Severe Acute Respiratory Syndrome)          |
| Vaccinia                                |                                                                          |
| Herpes Simplex Type I                   |                                                                          |
| Adenovirus Type 7 - at 8 oz. per gallon |                                                                          |

Enviro Care Neutral Disinfectant is effective against the following fungi according to the AOAC Fungicidal Test, modified in the presence of 5% horse serum (biostress load) and 400 ppm water hardness (calculated as CaCO3) at 20° C temperature. Dilution 1:64 (2 oz. per gallon of water).

- Aspergillus niger
- Trichophyton mentagrophytes
- Candida albicans

Enviro Care Neutral Disinfectant is mildewstatic as tested by CSMA, Fungal/Mold/Mildew Activity, Method #24. Modified in the presence of 5% fetal bovine biostress load and 400 ppm water hardness (calculated CaCO3) at 20° C temperature. Dilution 1:64 (2 oz. per gallon of water).

**For more information about RMC call 1.800.388.4762 or go to [www.rochestermidland.com](http://www.rochestermidland.com)**

